

## Explore the Possibilities

With the potential to generate a billion bases of DNA sequence per run, the Solexa Genome Analysis System provides researchers with the opportunity to sequence mammalian genomes in a matter of weeks rather than years, as is the case with present technologies. Furthermore, since our Clonal Single Molecule Array technology does not rely on expensive devices to achieve its unparalleled data density, the Solexa Sequencing-by-Synthesis approach is highly economical. Leveraging this technology, researchers can potentially resequence genomes for less than 1% of their current costs. With the capacity to accommodate up to 8 samples per flow cell, the Solexa Genome Analysis System can be tailored to the demands of your application. Whether your genome of interest is large or small, the Solexa Genome Analysis System is the ideal tool for your sequencing project.

## About Solexa

Solexa, Inc. (NASDAQ: SLXA) is developing and commercializing a new genetic analysis instrument system based on our novel Clonal Single Molecule Array™ technology and proprietary reversible terminator chemistry. This new platform, the Solexa Genome Analysis System, is being designed to offer a powerful new approach to some of today's most important genetic analysis applications, including:

- Whole Genome and Candidate Region resequencing
- Genome Wide Expression profiling
- Small RNA identification and quantification
- Hybridization array data validation

To learn more about Solexa, our technology, and our products, please visit [www.solexa.com](http://www.solexa.com) or email us at [info@solexa.com](mailto:info@solexa.com)

For more information on our Genomic Services, please contact us at:

**Solexa Inc.**  
Corporate Headquarters  
25861 Industrial Boulevard  
Hayward, CA 94545  
U.S.A.  
P: 510-670-9300  
F: 510-670-9302  
[info@solexa.com](mailto:info@solexa.com)  
[www.solexa.com](http://www.solexa.com)

**Solexa Ltd.**  
European Headquarters  
Chesterford Research Park  
Little Chesterford  
Essex, CB10 1XL  
United Kingdom  
P: +44 (0) 1799 532300  
F: +44 (0) 1799 532301

The Solexa Genome Analysis System leverages the company's proprietary Clonal Single Molecule Array™ technology and reversible terminator chemistry for generating rapid, economical, and accurate DNA sequence data. With the potential to yield up to 1 billion bases of DNA sequence per run at less than 1% of the cost of capillary-based methods, the Solexa Genome Analysis System is designed to enable researchers to dramatically improve the efficiency and speed of current applications and undertake new research that has been limited by lack of suitable technologies.

## THE WORKFLOW

- Random fragmentation of genomic DNA
- Addition of short adapters to the DNA fragments
- Immobilization of modified fragments to the flow cell
- Solid phase amplification to generate millions of distinct DNA clusters
- Base-by-base sequencing-by-synthesis using fluorescent reversible terminators
- Base calling and QC
- Application specific analysis

## Clonal Single Molecule Array™ Technology

Sequencing templates are immobilized on a proprietary flow cell surface designed to present the DNA in a manner that facilitates access to enzymes while ensuring high stability of surface-bound template and low non-specific binding of fluorescently labeled nucleotides. Solid phase amplification is employed to create up to 1,000 identical copies of each single molecule in close proximity (diameter of 1 micron or less). Since the process does not involve photolithography, mechanical spotting or positioning of beads into wells, the Clonal Single Molecule Array technology can achieve densities of up to 10 million single molecule clusters per square centimeter.

## Sequencing-by-Synthesis

Solexa's Sequencing-By-Synthesis (SBS) utilizes four proprietary fluorescently labeled modified nucleotides to sequence the millions of DNA clusters present on the flow cell surface. These nucleotides, specially designed to possess a reversible termination property, allow each cycle of the sequencing reaction to occur simultaneously in the presence of all four nucleotides (A, C, T, G). In each cycle, the polymerase is able to select the correct base to incorporate, with the natural competition between all four alternatives leading to higher accuracy than methods where only one nucleotide is present in the reaction mix at a time. Sequences where a particular base is repeated one after another (e.g., homopolymers) are addressed like any other sequence and with high accuracy.

## Analysis Pipeline

The Solexa sequencing approach is built around a very large number of short sequence reads. Deep sampling of more than ten-fold even coverage is required to generate a consensus and thus ensure high confidence in determination of genetic differences. Such differences are identified by comparison of sequence reads to a reference. Deep sampling allows the use of weighted "majority voting" and statistical analysis, similar to conventional methods, to identify homozygotes and heterozygotes and to distinguish sequencing errors. Each raw read base has an assigned quality score so that the software can apply a weighting factor in calling differences and generating confidence scores.

The suite of software from Solexa will enable users to align sequences to a reference in resequencing applications. Developed in collaboration with leading researchers, Solexa's software suite includes the full range of data collection, processing, and analysis modules to streamline collection and analysis of data with minimal user intervention. The open format of the software allows for easy access to the data at various stages of processing and analysis using simple application program interfaces.

# Sequencing Technology Overview

